

STERILISED SPOONS



Where sterile techniques are required for collection of samples, these sterilised spoons are ideal for use with Whirl-Pak® Sterilised Closure bags.

These easy-to-use spoons are made of sturdy, durable polystyrene.

Instructions on how to use

The procedure for contamination-free use is to open the bag, push one spoon up, pull it out, and re-close the bag.

Spoons are packed 20 per bag. Box of 100.



PRODUCT CODE	DESCRIPTION	PACK SIZE
B01041-A-WA	Spoon Disposable 15cm - Sterilised	PK100

CELLULOSE SPECI-SPONGE® BAGS

Whirl-Pak® Cellulose Speci-Sponge® bags are designed for environmental surface sampling of work areas, equipment, animal carcasses, and any other place where testing for *Listeria*, *Salmonella*, *E.coli*, coliforms, and other foodborne pathogens, or anthrax is required.

Each sponge is packed in a sterile, sealed write-on bag with a dry sponge.



PRODUCT CODE	DESCRIPTION	PACK SIZE
B01245WA	Bag WhirlPak Speci-Sponge 115 x 230mm (WxL) - 532mL	PK100

SODIUM THIOSULPHATE STERILISED BAGS

An ideal sustainable shipping and storage solution, specifically designed for sampling chlorinated water.

Suitable in all sectors including municipalities, swimming pools, waste water and other areas to test for fecal coliforms, legionella and other bacteria. These Thio Sterilised bags provide a secure, contaminant free flexible container that ensures dependable analysis results.

Features:

- Maximum protection against leaks with their signature perimeter seal
- Write on strip for easy sample identification
- Extreme film clarity and strength
- Self standing bag
- RNase, DNase, Pyrogen and BPH free
- Aseptic design with Tyvek® opening pull tabs - this feature keeps your hands away for a safe, aseptic collection technique each time
- Supplied sterile - the bags go through post-manufacturing sterilisation




Sodium Thiosulphate tablet included in bag

Learn about how they work



Store and transport samples while reducing costs, space and waste



Space saving - No need to find storage for bulky boxes!

*Containers and boxes are for illustration purposes only. Not to scale.

PRODUCT CODE	DESCRIPTION	PACK SIZE
B01040WA	Bag WhirlPak Sodium Thiosulphate *Standard Bag 75 x 185mm (WxL) - 100mL <small>CLEARANCE</small>	PK100
B01402WA	Bag WhirlPak Sodium Thiosulphate Stand Up Bag 75 x 185mm (WxL) - 100mL	PK100
B01403WA	Bag Whirlpak Sodium Thiosulphate Stand Up Bag 115 x 230mm (WxL) - 300mL <small>CLEARANCE</small>	PK100