

OVERSTOCK SALE

Hurry! Stock up and save big on these items until the end of April, or while stock lasts!

METER AQUA CPA (pH/COND/TEMP)

Designed and built in Australia, the Aqua-Series is used from mines to meat testing, councils to construction and water treatment to wine making.

- pH
- Conductivity (TDS and Salinity)
- Temperature
- Sensors and buffers included
- Ready to use out of the box
- Splash proof
- 9V battery
- Fold out stand
- Australian Made



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
121135TPS-EC1001	Meter Aqua-CPA Cond/TDS/Salinity/pH/Temp 1m Cable	Each	Login	Login

Available Stock: 3

ELECTRODE STAND/HOLDER

Flexible arm type stand for holding electrodes

- Holds 2 x 12mm and 1 x 6.5mm diameter electrodes
- Heavy base for stability




PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
130088TPS	Stand for Electrode TPS - Flexible arm type	Each	Login	Login

Available Stock: 2

TUBE CENTRIFUGE 50ML, STERILE & RACKED

Made from polypropylene these tubes are certified non-pyrogenic and DNase/RNase Free

- Supplied in cardboard racks (25)
- Conical bottom (29mm x 115mm long)
- Autoclavable to 122°C
- Graduated from 5mL to 45mL x 2.5mL
- Max RCF 9,500g
- Flat cap
- Sterile in inner packs of 25 (racked)



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
650.550.250	Tube Centrifuge 50mL PP Racked	PK500	Login	Login

*Bulky charge may apply

Available Stock: 120

TROLLEY ADJUSTABLE HEIGHT 2 SHELVES

- Overall dimensions 900L x 500W x 970H mm
- 100mm casters (two with brakes)
- High strength plastic shelves 715 x 455mm
- Easy to assemble
- Load capacity: 275kg



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
800.011.020	Trolley Adjustable Height 2 Shelves	Each	Login	Login

*Bulky charge may apply

Available Stock: 12


OVERSTOCK SALE



Swabs


APPLICATOR

- Wooden shaft
- No tip
- Non sterile
- 150mm x 2mm




WOODEN SHAFT

- Wooden shaft
- Rayon tip
- Non sterile
- 150mm x 2mm



PLASTIC SHAFT

- Plastic shaft (PP)
- Rayon tip
- Non sterile
- 150mm x 2mm



PRODUCT CODE	PACK SIZE	PRICE	SALE
710.050.005	PK1000	Login	Login

Available Stock: 150

PRODUCT CODE	PACK SIZE	PRICE	SALE
710.050.020	PK100	Login	Login

Available Stock: 529

PRODUCT CODE	PACK SIZE	PRICE	SALE
710.050.021	PK100	Login	Login

Available Stock: 68

TRANSPORT PLASTIC

- Swab in tube
- Rayon tip
- Sterile peel pack, individually wrapped
- PP tube 13mm x 150mm



TRANSPORT AMIES MEDIA

- Amies (without charcoal) media
- Swab in tube
- Rayon tip
- Sterile peel pack, individually wrapped
- PP Tube 13mm x 150mm



**SHORT EXPIRY DATE
EXPIRES JULY 2024**

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE	QTY AVAIL
710.050.061	Swab transport plastic (PP) shaft sterile	PK100	Login	Login	166
710.050.120	Swab transport plastic shaft amies media *Short expiry date, expires July 2024*	PK100	Login	Login	72

OVERSTOCK SALE

SLEEVE EXTENDERS

- Disposable
- Polyethylene
- 40cm
- 3GSM



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
900.200.005	Sleeve Extenders Disposable	PK100	Login	Login

Available Stock: 174

BEARD COVERS

- Disposable
- White
- 12GSM
- One loop (behind head)



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
900.150.005	Beard Cover Disposable White	PK100	Login	Login

Available Stock: 236

STOPPERS CONICAL PP

- Made from polypropylene



Multi-Buy minimum 10 pcs

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE	QTY AVAIL
650.015.210	Stopper Conical PP NS 10/19 - Polypropylene	Each	Login	Login	317
650.015.212	Stopper Conical PP NS 12/21 - Polypropylene	Each	Login	Login	171
650.015.214	Stopper Conical PP NS 14/23 - Polypropylene	Each	Login	Login	61
650.015.219	Stopper Conical PP NS 19/26 - Polypropylene	Each	Login	Login	114
650.015.224	Stopper Conical PP NS 24/29 - Polypropylene	Each	Login	Login	135
650.015.229	Stopper Conical PP NS 29/32 - Polypropylene	Each	Login	Login	51

Prices exclude GST. Promo validity date (30/04/2024). E.&O.E. Stock availability is current as of 9/4/24.

OVERSTOCK SALE



RACK CARDBOARD FREEZER

- All feature alphanumeric grid on lid for tube ID



650.100.100

- For tubes: 15mL (17mm diameter)
- 134 x 134 x 115 mm (LxWxH)
- 7 x 7 configuration



650.100.102

- For tubes: 30mL (30mm diameter)
- 134 x 134 x 122 mm (LxWxH)
- 4 x 4 configuration



650.100.104

- For tubes: 15mm diameter
- 134 x 134 x 47 mm (LxWxH)
- 8 x 8 configuration



650.100.108

- For tubes: 11mm diameter
- 134 x 134 x 47 mm (LxWxH)
- 10 x 10 configuration

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE	QTY AVAIL
650.100.100	Rack Cardboard Freezer 49 Place (7x7)	Each	Login	Login	221
650.100.102	Rack Cardboard Freezer 16 Place (4x4)	Each	Login	Login	221
650.100.104	Rack Cardboard Freezer 64 Place 47mm H	Each	Login	Login	471
650.100.108	Rack Cardboard Freezer 100 Place 47mm H	Each	Login	Login	1,163

*Bulky charge may apply

SLIDE BOX 100 PLACE

- For slides 25 x 75mm
- Made from high impact ABS
- Cork lined
- Numbered slots & slide inventory sheet for slide ID
- 199 x 170 x 33 mm (LxWxH)



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
650.100.032	Slide Box 100 Place Blue	Each	Login	Login

Available Stock: 1,013

SLIDE MAILER CARDBOARD 20 PLACE

- For slides 25 x 75mm
- Made from heavy weight cardboard
- Numbered compartments & index sheet for slide ID
- 334 x 205 mm (LxW)



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
650.100.066	Slide Mailer Cardboard 20 Place	PK5	Login	Login

Available Stock: 172

OVERSTOCK SALE

RACK 80 WELL

- For tubes: 1.5mL & 2mL
- Made from heavy duty polypropylene
- Temperature range: -80 to 121°C
- Moulded numbered on the surface for tube ID
- Autoclavable 121°C
- 228 x 66 x 22 mm (LxWxH)
- 5 x 16 configuration



650.100.254



650.100.258

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
650.100.254	Rack 80 Well Pink	Each	Login	Login
Available Stock: 48				
650.100.258	Rack 80 Well Orange	Each	Login	Login
Available Stock: 32				

RACK 96 WELL

- For tubes: 1.5/2mL & 0.5/07mL
- Made from heavy duty polypropylene
- Temperature range: -80 to 121°C
- Moulded numbered on the surface for tube ID
- Autoclavable up to 121°C
- 213 x 111 x 25 mm (LxWxH)
- 8 x 12 Configuration



650.100.260




650.100.268

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
650.100.260	Reversible Well Orange	Each	Login	Login
Available Stock: 97				
650.100.268	Reversible Well Blue	Each	Login	Login
Available Stock: 104				

RACK 100 PLACE

- For tubes: 1.5mL & 2mL
- Made from polypropylene
- Temperature range: -80 to 121°C
- Alphanumeric grid on hinge lid for tube ID
- 130 x 130 x 45 mm (LxWxH)
- 10 x 10 configuration



650.100.222 650.100.224

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
650.100.222	Rack 100 Place Blue Tube 1.5mL - 2mL	Each	Login	Login
Available Stock: 60				
650.100.224	Rack 100 Place Green Tube 1.5mL - 2mL	Each	Login	Login
Available Stock: 32				

RACK 4-WAY YELLOW

- For tubes: 50mL x 4, 15mL x 12, 1.5mL x 32, 0.5mL x 32
- Made from heavy duty polypropylene
- Autoclavable up to 121°C
- Inter-lockable design to customise your storage requirements
- 174 x 95 x 52 mm (LxWxH)



**YELLOW
RACKS ON
CLEARANCE!**



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
650.100.208	Rack 4-Way Yellow	Each	Login	Login

Available Stock: 49

RACK SNAP TOGETHER 50 PLACE TUBE 15/50ML

- For tubes: 15mL (17mm D) x 30 and 50mL (28mm D) x 20
- Easy to assemble (supplied flat)
- 207 x 170 x 60 mm (LxWxH)
- 6 x 5 config (50mL holes), 5 x 5 (15mL holes)
- Designed for stable use in water baths



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
650.100.121	Rack Snap Together 50 Place Tube 15/50mL	Each	Login	Login

Available Stock: 855