

Incubator



INC 125 F DIGITAL

NEW PRODUCT

The INC 125 F digital convection incubator is ideal for microbiological applications up to 100°C with precise temperature control. The supply of preheated air into the temperature chamber ensures fast, stable and uniform temperature control.

Features + Benefits

- Easy to clean
- Optimal temperature control
- Simple and intuitive operation
- Versatile accessories

Specifications:

- 125L capacity
- Temperature range: RT +8°C to 100°C
- Temperature stability: ±0.3K
- Dimensions: 700 x 720 x 635mm
- 2x chrome-plated wire grid trays included
 - can fit up to 6 grid shelves



Inner glass door for easy observation without temperature loss

See the IKA Incubator in action here:



WATCH NOW

Large LCD display for easy reading



Counter, timer and auto-timer function



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
00201171721	Incubator IKA INC 125 F Digital	Each	LOGIN	LOGIN

Coveralls



Ansell BioClean™ Coveralls are designed for cleanroom and controlled environments, offering superior protection and comfort.

Made from low-linting, antistatic material, they provide an effective barrier against particles and liquid splashes. Lightweight and breathable, these coveralls help ensure contamination control while maintaining comfort.



Compliance & Quality

BioClean™ products offer innovative solutions which are designed to the highest quality that meet user needs.



Regulatory Standards

BioClean™ products are rigorously tested and fully comply with industry standards and regulations.



Cleanroom Solutions

High quality cleanroom compatible solutions processed & packed in ISO Class 4 & 5 cleanrooms.

STERILE

Hood and Integrated Booties

Optimum Protection



Features + Benefits

Sterile, disposable coverall with hood and integrated boots

Reduced Contaminant Risk

- Lightweight CleanTough material means the coverall is low-linting and ISO cleanroom compliant, reducing contamination risks

Optimised Protection

- Front zip with a protective flap, resisting liquid splashes and chemicals
- Thumb loops and an elasticated hood, back, cuffs and ankles ensure a secure fit

Enhanced Features

- This antistatic coverall features integrated boots and slip-resistant soles, for added workplace safety in laboratory settings

Sterility Assurance Level

- Special packaging eliminates cross contamination, while gamma irradiation sterilisation ensures a minimum SAL of 10^{-6}



Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Pharmaceutical Manufacturing
- Biotechnology Manufacturing
- Medical Device Manufacturing

PRODUCT CODE	DESCRIPTION	SIZE	PACK SIZE	PRICE
S-BDFC-XS	Coverall with Hood and Integrated Booties	XS	PK20	LOGIN
S-BDFC-S	Coverall with Hood and Integrated Booties	S	PK20	LOGIN
S-BDFC-M	Coverall with Hood and Integrated Booties	M	PK20	LOGIN
S-BDFC-L	Coverall with Hood and Integrated Booties	L	PK20	LOGIN
S-BDFC-XL	Coverall with Hood and Integrated Booties	XL	PK20	LOGIN
S-BDFC-XXL	Coverall with Hood and Integrated Booties	2XL	PK20	LOGIN

Sizing guide found on individual product page

Coveralls

STERILE



Hood

Features + Benefits

Sterile, disposable coverall with hood

Reduced contaminant risk

- Featuring a hood and front zip with sealable cover and a protective flap, reducing the risk of contaminant leakage through the seams
- Manufactured in compliance with ISO cleanroom standards

Specialised protection

- Manufactured from CleanTough material to deliver a low-linting, antistatic, sterile coverall that protects against a range of chemicals

Optimised fit

- Equipped with an elasticated hood, back, cuffs and ankles, as well as thumb loops to ensure a secure hold and enhanced fit

Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Pharmaceutical Manufacturing
- Biotechnology Manufacturing
- Medical Device Manufacturing



*Sterile shoe covers and gloves sold separately

PRODUCT CODE	DESCRIPTION	SIZE	PACK SIZE	PRICE
S-BDCHT-XS	Coverall BioClean with Hood - Sterile	XS	PK20	LOGIN
S-BDCHT-S	Coverall BioClean with Hood - Sterile	S	PK20	LOGIN
S-BDCHT-M	Coverall BioClean with Hood - Sterile	M	PK20	LOGIN
S-BDCHT-L	Coverall BioClean with Hood - Sterile	L	PK20	LOGIN
S-BDCHT-XL	Coverall BioClean with Hood - Sterile (Expiry Date 01/05/25)	XL	PK20	LOGIN
S-BDCHT-2XL	Coverall BioClean with Hood - Sterile	2XL	PK20	LOGIN
S-BDCHT-3XL	Coverall BioClean with Hood - Sterile	3XL	PK15	LOGIN
S-BDCHT-4XL	Coverall BioClean with Hood - Sterile	4XL	PK15	LOGIN

Sale Price
Close Expiry

Sizing guide found on individual product page

Need Sterile Shoe Covers?

[BUY HERE](#)

Coverall BioClean with Hood - Sterile
Size XL Expiry Date 01/05/2025

[BUY NOW](#)

Coveralls

NON-STERILE



Hood

Features + Benefits

Non-sterile, disposable coverall with hood

Reduced contaminant risk

- Featuring a hood and front zip with sealable cover and a protective flap, reducing the risk of contaminant leakage through the seams
- Manufactured in compliance with ISO cleanroom standards

Specialised protection

- Manufactured from CleanTough material, granting the garment lightweight and low-linting characteristics

Optimised fit

- Equipped with an elasticated hood, back, cuffs and ankles, as well as thumb loops to ensure a secure hold and enhanced fit

Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Pharmaceutical Manufacturing
- Biotechnology Manufacturing
- Medical Device Manufacturing



*Shoe covers and gloves sold separately

PRODUCT CODE	DESCRIPTION	SIZE	PACK SIZE	PRICE
BDCHT-XS	Coverall BioClean with Hood - Non-Sterile	XS	PK20	LOGIN
BDCHT-S	Coverall BioClean with Hood - Non-Sterile	S	PK20	LOGIN
BDCHT-M	Coverall BioClean with Hood - Non-Sterile	M	PK20	LOGIN
BDCHT-L	Coverall BioClean with Hood - Non-Sterile	L	PK20	LOGIN
BDCHT-XL	Coverall BioClean with Hood - Non-Sterile	XL	PK20	LOGIN
BDCHT-2XL	Coverall BioClean with Hood - Non-Sterile	2XL	PK15	LOGIN
BDCHT-3XL	Coverall BioClean with Hood - Non-Sterile	3XL	PK15	LOGIN
BDCHT-4XL	Coverall BioClean with Hood - Non-Sterile	4XL	PK15	LOGIN

Sale Price
Clearance

Sizing guide found on individual product page

Need Non-Sterile Shoe Covers?

[BUY HERE](#)