

# Homogeniser

## T25 HOMOGENISER PACKAGE

High performance dispersing instrument for volumes from 1–2,000mL (H<sub>2</sub>O) with digital speed display. It offers a wide speed range from 3,000–25,000rpm that enables users to work at high circumferential speeds even with small rotor diameters.

### Features + Benefits

- Digital speed display
- Electronic speed control
- Electronic overload protection
- Stainless steel dispersing elements can be cleaned quickly and easily
- Error code display

### Supplied with:

- T25 Digital
- S25N/18G dispersing tool
- R182 bosshead clamp
- R1826 stand



**lab**  
co  
SCIENTIFIC

### Applications

- Food & Beverage
- Pharmaceutical/Biotechnology
- Cosmetics & Personal Care
- Dairy
- Sample Preparation
- Scientific Research
- Environmental Analysis



| PRODUCT CODE                | DESCRIPTION             | PACK SIZE | PRICE                 | SALE                  |
|-----------------------------|-------------------------|-----------|-----------------------|-----------------------|
| <a href="#">00090242001</a> | Homogeniser T25 Package | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |

# Overhead Stirrers

## LABCO 20L OVERHEAD STIRRER

This overhead stirrer can facilitate the transition from small scale experimentation to large scale production.

### Features + Benefits

- Up to 20L stirring capacity
- Push-through stirring shaft
- Overload protection
- Wide, accurate speed range
- Smooth start avoids splashing
- DC brushless motor
- LCD displays set and actual speed/torque



\*20L Overhead Stirrer pictured with stand, rod, bosshead & 4 blade paddle. Beaker sold separately.

| PRODUCT CODE                | DESCRIPTION                        | PACK SIZE | PRICE                 | SALE                      |
|-----------------------------|------------------------------------|-----------|-----------------------|---------------------------|
| <a href="#">400.103.020</a> | Overhead Stirrer 20L               | Each      | <a href="#">LOGIN</a> | -                         |
| <a href="#">400.103.010</a> | Stand including Bosshead Clamp     | Each      | <a href="#">LOGIN</a> | <b>FREE WITH PURCHASE</b> |
| <a href="#">400.030.050</a> | Paddle 4 Blade 316 Stainless Steel | Each      | <a href="#">LOGIN</a> | <b>FREE WITH PURCHASE</b> |

### FREE WITH PURCHASE OF STIRRER

1x Stand, including bosshead clamp  
1x 4 Blade paddle



# Overhead Stirrers

Powerful, mechanically controlled laboratory stirrers designed for highly viscous applications. The stirrers of the RW series are suitable for intensive mixing for use in laboratories and pilot plants.

## RW20 OVERHEAD STIRRER PACKAGE

### Features + Benefits

- The bestseller in the laboratory
- Digital display
- Robust, slimline, ergonomic design
- Constant power-drive
- Two speed ranges for universal use
- Push-through agitator shafts (only when stationary)

### Technical Specifications:

- Max. Stirring Volume: 20L (H<sub>2</sub>O)
- Max. Viscosity: 10,000mPas
- Speed Range: 60 – 2,000rpm
- Max. Torque: 150Ncm

### Supplied with:

- RW20 Stirrer
- R182 Bosshead Clamp
- R1826 Stand Plate
- R1342 Paddle

| PRODUCT CODE                | DESCRIPTION  | PACK SIZE | PRICE                 | SALE                  |
|-----------------------------|--------------|-----------|-----------------------|-----------------------|
| <a href="#">0009011600I</a> | RW20 Package | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |



\*Beaker clamp not included with delivery.

## RW28 OVERHEAD STIRRER PACKAGE

### Features + Benefits

- Digital speed display
- Infinitely adjustable speed
- Push-through agitator shafts
- Overload protection
- Error code display
- Robust, ergonomic design
- Quiet operation
- Constant power-drive

### Technical Specifications:

- Max. Stirring Volume: 80L (H<sub>2</sub>O)
- Max. Viscosity: 50,000mPas
- Speed Range: 60 – 1,400rpm
- Max. Torque: 900Ncm

### Supplied with:

- RW28 Stirrer
- R271 Bosshead Clamp
- R2723 Telescopic Stand
- R1385 Paddle

| PRODUCT CODE                | DESCRIPTION  | PACK SIZE | PRICE                 | SALE                  |
|-----------------------------|--------------|-----------|-----------------------|-----------------------|
| <a href="#">0009023200I</a> | RW28 Package | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |

Price valid for  
stock only



# Mills

## BASIC MILL A11

### Features + Benefits

- Ideal for impact grinding and cutting grinding
  - 80mL capacity
  - 28,000rpm fixed speed
  - Moist and gluey materials can be pulverised by adding water
  - Grinding chamber made of Tefcel (PTFE, glass fiber-reinforced) with stainless steel inlet (1.4571)
  - Overload protection
  - 1.5kg weight
  - A11.1 beater\* included with delivery
- \*Not suitable for wet milling



## BASIC MILL A10

### Features + Benefits

- Ideal for grinding hard, brittle, soft and fibrous materials
- 50mL capacity
- 25,000rpm variable speed, deviation  $\pm 500$ rpm
- Digital timer
- Counter (Displays grinding time)
- Integrated cooling connectors
- Interval function
- Brushless motor
- Error code display
- 3.13kg weight
- A10.1 cutter included with delivery



### Applications

- Agriculture Industry - grain, seeds
- Brewery - hop, malt pellets
- Building materials
- Detergents, mega pearls
- Energy Industry - coal
- Animal feed
- Food industry - rice, spices, bread, meat
- Medical Industry - tablets, pills
- Paint Industry - colour pigments, resins
- Paper Industry - Paper, filter, cellulose
- Pharmacy - Vitamin tablets, tea leaves



| PRODUCT CODE               | DESCRIPTION        | PACK SIZE | PRICE                 | SALE                  |
|----------------------------|--------------------|-----------|-----------------------|-----------------------|
| <a href="#">0002900000</a> | A11 IKA Basic Mill | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0004020700</a> | A10 IKA Basic Mill | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |

# Viscometers

The ROTAVISC series determines the liquid viscosity in all areas of application ranging from the laboratory to quality control.

Versatile and precise viscometers are essential for industries where viscosity measurements are critical, such as pharmaceuticals, chemicals, food production and cosmetics.



**IKA®**

## Features + Benefits

- Parameter display
- Precise results From 0.01 – 200rpm
- Stepless measurement
- Stepless positioning
- Temperature control
- Three measuring modes

## Supplied with:

- Viscometer
- Rotastand
- Relevant spindle set
- PT100.8 temperature sensor
- Guard rail
- USB cable



## APPLICATIONS


### Personal Care/Cosmetics

- Hair Care Products
- Fragrance Productions
- Toothpaste Manufacturing
- Shampoo & Conditioner
- Sunscreen & Skincare Products



### Pharmaceutical

- Pharmaceutical Formulations
- Oral Care Products
- Gels for Wound Care
- Parenteral Formulations
- Quality Control



### Food & Beverage

- Sauces and Gravies
- Ketchup and Mustard
- Chocolate
- Batters
- Dairy Products
- Smoothies
- Syrups



### Paints, Coatings, Lubricants, Paper & More

- Paint Formulation
- Coating Thickness
- Adhesive Production
- Grease Manufacturing
- Pulp Viscosity
- Ink Manufacturing



### Unsure which model will suit your applications?

Simply complete [this form](#) and our experts will be in touch.

| PRODUCT CODE                | DESCRIPTION       | VISCOSITY RANGE        | PACK SIZE | PRICE                 |
|-----------------------------|-------------------|------------------------|-----------|-----------------------|
| <a href="#">0025000310I</a> | ROTAVISC lo-vi    | 1 to 6,000,000mPas     | Each      | <a href="#">LOGIN</a> |
| <a href="#">0025000311I</a> | ROTAVISC me-vi    | 100 to 40,000,000mPas  | Each      | <a href="#">LOGIN</a> |
| <a href="#">0025000312I</a> | ROTAVISC hi-vi I  | 200 to 80,000,000mPas  | Each      | <a href="#">LOGIN</a> |
| <a href="#">0025000313I</a> | ROTAVISC hi-vi II | 800 to 320,000,000mPas | Each      | <a href="#">LOGIN</a> |



# Magnetic Hotplate Stirrer



Experience the convenience of simultaneously heating and stirring solutions in your lab. Resistant to corrosion, this unit is easy to clean and has a digital LCD display that shows both set & actual temperature and speed.

### Features + Benefits

- 20L stirring capacity
- Speed range 100-1,500rpm
- Three selectable heating modes: fast, standard & stable
- HOT warning, indicating residual heat
- Digital LCD displays both set & actual temperature & speed
- Maintenance free DC Brushless motor
- Automatic recognition of PT1000 once connected
- External PT1000 temperature sensor improves accuracy to  $\pm 0.2^{\circ}\text{C}$  degrees



### Supplied with:

- Protective cover (400.100.208)
- Power cable
- User manual



Stirrer pictured with protective cover

| PRODUCT CODE                | DESCRIPTION                     | PACK SIZE | PRICE                 | SALE                  |
|-----------------------------|---------------------------------|-----------|-----------------------|-----------------------|
| <a href="#">400.100.112</a> | Magnetic Hotplate Stirrer 20L   | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">400.100.218</a> | Sensor Temperature PT1000 230mm | Each      | <a href="#">LOGIN</a> | FREE WITH PURCHASE    |
| <a href="#">400.100.225</a> | Clamp Support for PT1000 Sensor | Each      | <a href="#">LOGIN</a> | FREE WITH PURCHASE    |

# Pipettors

## IKA® PETTE VARIO

An ergonomic feel combined with uncompromising precise performance.  
The IKA® pipettes can be used for many applications in the modern laboratory.

### Ergonomic design

Easy one-handed operation with exchangeable grips for a comfortable fit

### Colour coding

Simplifies the rapid selection of a suitable pipette

### Durable construction

Made with high quality materials, resistance to impact, UV and chemicals, ensuring long-term reliable dosing

### Thermally insulated

Reduces volume errors

### Fully autoclavable

Ensuring easy sterilisation. This guarantees fault-free, precise and safe operation

### User-friendly calibration

Easy-to-use IKA® software for planning, performing and recording calibration tests

### Simple maintenance


The IKA® Pettes can be disassembled, cleaned and adjusted in just a few simple steps without any additional tool

### Made in Germany

100% produced, assembled and calibrated in Germany



Supplied with 3 ergonomic grips



### Pipettor | Single Channel | IKA®



| PRODUCT CODE                | VOLUME RANGE  | PACK SIZE | PRICE                 | SALE                  |
|-----------------------------|---------------|-----------|-----------------------|-----------------------|
| <a href="#">0020011210I</a> | 0.1 - 2µL     | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0020011211I</a> | 0.5 - 10µL    | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0020011213I</a> | 2 - 20µL      | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0020011214I</a> | 10 - 100µL    | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0020011215I</a> | 20 - 200µL    | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0020011216I</a> | 100 - 1,000µL | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0020011217I</a> | 0.5 - 5mL     | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0020011218I</a> | 1 - 10mL      | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |

### Accessories

| PRODUCT CODE                | DESCRIPTION                | PACK SIZE | PRICE                 | SALE                  |
|-----------------------------|----------------------------|-----------|-----------------------|-----------------------|
| <a href="#">0020018987I</a> | Pipette Rack Linea         | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0020107134I</a> | Pipette Stand Ronda        | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0020016449I</a> | Filter for Pette 0.5 - 5mL | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |
| <a href="#">0020016450I</a> | Filter for Pette 1 - 10mL  | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |



Ronda Stand  
0020107134I



Linea Stand  
0020018987I

# pH Meters



## ZENTEST® SMART | POCKET TESTERS

This Apera 60-Z Smart Pocket Test is the next-gen tester designed for pH measurements of food/solid/viscous samplings such as cheese, meat, fruit, dough, yogurt, sushi rice, jam, sauce as well as skincare products.

Powered by the ZenTest® Mobile App and connected via Bluetooth®. Save and share your measurement and calibration data right on the palm of your hand via email. Data will never get lost thanks to the cloud-based data system.

### Features + Benefits

- Spear pH probe designed for food and solid samples
- Smart self diagnosis
- Upgraded sensor technology
- Replaceable electrode
- Durable design IP67 waterproof and dustproof
- Cloud based data management
- Hybrid design (can be used with or without a smart phone)

| PRODUCT CODE |             | PH60S-Z SPEAR   |
|--------------|-------------|---|
| PARAMETER    |             |   |
| pH           | RANGE       | -2.00 to 16.00 pH   |
|              | RESOLUTION  | ±0.01pH; ±1 digit   |
|              | CALIBRATION | 1 to 3 points Auto. calibration (recognises 5 types of standards) |
| ORP (mV)     | RANGE       | ±1,000 mV*  |

\*optional with [ORP60-DA](#) probe

| PRODUCT CODE            | DESCRIPTION                  | PACK SIZE | PRICE                 | SALE                  |
|-------------------------|------------------------------|-----------|-----------------------|-----------------------|
| <a href="#">PH60S-Z</a> | ZenTest® Spear Pocket Tester | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |

## 850 VALUE SERIES | HANDHELD

### Low cost, reliable, portable meter with simple operation

- Large LCD display made from acrylic and is resistant to scratches
- Quick and easy calibration
- IP57 waterproof rating
- Comes with everything you need, including calibration solutions and electrodes, in a convenient carry case

| PRODUCT CODE          | DESCRIPTION | PACK SIZE | PRICE                 | SALE                  |
|-----------------------|-------------|-----------|-----------------------|-----------------------|
| <a href="#">PH850</a> | pH Meter    | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |



# pH Meters

## 910 LAB GRADE SERIES | BENCHTOP

This 910 Series pH benchtop meter is designed for laboratory grade measurement for general water solutions.

It features an intuitive operating system, automatic calibration with smart guidance and self-diagnosis, automatic data logger, and many more features while retaining simplicity in use.

- TFT colour screen with graphics and text guides
- Supplied with flexible electrode holder
- Complete kit with ready-to-use buffers
- IP54 waterproof rating
- Stable reading indication with a smiley face



| PRODUCT CODE          | DESCRIPTION       | PACK SIZE | PRICE                 | SALE                  |
|-----------------------|-------------------|-----------|-----------------------|-----------------------|
| <a href="#">PH910</a> | pH Benchtop Meter | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |

## 9500 RESEARCH GRADE SERIES | BENCHTOP

This 9500 Series pH benchtop meter is designed for research laboratory professionals, offering long term reliability and accuracy in measurement.

It is supplied with a 'TestBench' that includes a flexible electrode holder, an intelligent magnetic stirrer, and a buffer organiser.

- GLP data logging, data export (to PC or Printer)
- Calibration password protection
- Buffer organiser
- Self-diagnosis
- USB keyboard
- Supplied with flexible electrode holder
- Supplied with electrodes
- Includes magnetic stirrer



| PRODUCT CODE           | DESCRIPTION         | PACK SIZE | PRICE                 | SALE                  |
|------------------------|---------------------|-----------|-----------------------|-----------------------|
| <a href="#">PH9500</a> | pH Test Bench Meter | Each      | <a href="#">LOGIN</a> | <a href="#">LOGIN</a> |

## INTERMEDIATE JUNCTION pH ELECTRODE

- Renewable intermediate junction
- Premium performance
- Super long life
- Heavy duty

### Specifications:

- 0 to 14pH
- Glass spear membrane
- 0 to 60°C temperature range

### Applications:

Food, beverage and penetration



| PRODUCT CODE                 | DESCRIPTION                        | PACK SIZE | PRICE                 |
|------------------------------|------------------------------------|-----------|-----------------------|
| <a href="#">IJ44A-BNC-1M</a> | Electrode pH Intermediate Junction | Each      | <a href="#">LOGIN</a> |