

IKA® BLITZ IS BACK!

SAVE ON EQUIPMENT

- Overhead Stirrers
- Magnetic Stirrer Hotplates
- Pipettors
- Ovens and Incubators
- Homogenisers

LIMITED TIME OFFER

CLEARANCE SPECIALS

- Save over \$2,000 on IKA Mill 100 Control!
- IKA I-Mag + SP300 Frame Package - Save BIG!
- Savings of up to 46% inside!



Overhead Stirrers

RW20 OVERHEAD STIRRER

Robust, slimline, ergonomic design

- Digital speed display
- Robust, slimline, ergonomic design
- Max. stirring volume: 20L (H₂O)
- Max. viscosity: 10,000mPas
- Speed range: 60 – 2,000rpm
- Max. torque: 150Ncm
- Constant power-drive
- Push-through agitator shafts (when stationary)



Supplied with:

- RW20
- Check key

RW20 OVERHEAD STIRRER PACKAGE

Complete stirring setup

- Digital speed display
- Robust, slimline, ergonomic design
- Max. stirring volume: 20L (H₂O)
- Max. viscosity: 10,000mPas
- Speed range: 60 – 2,000rpm
- Max. torque: 150Ncm
- Constant power-drive
- Push-through agitator shafts (when stationary)

Supplied with:

- RW20
- R182 bosshead
- R1826 plate stand
- R1342 paddle



*RH3 Clamp supplied separately.

PRODUCT CODE	DESCRIPTION	PRICE	SALE
0003593000I	Overhead Stirrer RW20	LOGIN	LOGIN

PRODUCT CODE	DESCRIPTION	PRICE	SALE
0009011600I	RW20 Package	LOGIN	LOGIN

RW28 OVERHEAD STIRRER PACKAGE

Robust, ergonomic design for larger volumes

Powerful, mechanically controlled laboratory stirrer designed for highly viscous applications for quantities up to 80L (H₂O). It is suitable for intensive mixing for use in laboratories and pilot plants.

The stirrer comes equipped with a protection for overheating of the motor by means of self-locking temperature limiter. Additionally, to secure bowls, a stirring shaft protection and clamping fixture are optionally available.

- Digital speed display
- Quiet operation
- Max. stirring volume: 80L (H₂O)
- Max. viscosity: 50,000mPas
- Infinitely adjustable speed of 60 - 1,400rpm (in two speed ranges)
- Max. torque: 900Ncm
- Constant power-drive
- Push-through agitator shafts (when stationary)

Supplied with:

- RW28
- R271 bosshead
- R2723 telescopic stand
- R1385 paddle



IKA®

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0009023200I	RW28 Package	Each	LOGIN	LOGIN

Hotplate Stirrer

RCT BASIC HOTPLATE STIRRER

NEW!

The latest generation of RCT basic presents itself in proven quality and with numerous new features

- Tempered, shatterproof glass surface for optimal chemical resistance and easy cleaning
- Large, easy-to-read LED display
- Displays set and actual temperature, heating status and temperature sensor
- USB and RS232 interface
- Compatible with labworldsoft® laboratory software
- QR code for quick access to device information, accessories, downloads and extended warranty
- Easily accessible main switch on the front of the device

Proven technology

- Heating plate made of aluminum for optimal and homogeneous heat transfer
- Excellent magnetic coupling
- Soft start prevents magnetic stir bars from breaking off during the start phase
- Two optimized temperature control modes ensure fast heating or precise temperature control without overshooting
- Push and turn buttons for independent adjustment of the set points and starting / stopping of temperature and speed

Safe, robust and compliant with standards

- In accordance with the DIN EN 61010-2-010 standard

Specifications:

- Plate diameter: 135mm Ø
- Stirring volume: 20L (H₂O)
- Speed range: 50 - 1,500rpm
- Stirrer bar length: 20 - 80mm
- Heating temp. range: Ambient - 310°C
- Dimensions: 160 x 85 x 270mm
- Weight: 2.5kg

Supplied with:

- RCT Basic
- 30mm magnetic stirrer bar
- 40mm magnetic stirrer bar
- Screw driver
- PT1000.60 temperature sensor, stainless steel



*Thermometer, beaker, stand and clamp not included.

WATCH NOW

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
00250059271	RCT Basic Hotplate Stirrer	Each	LOGIN	LOGIN

SAVE
28%!

Pipettors

UPGRADE YOUR PIPETTING

An ergonomic design, combined with uncompromising precise performance. The IKA® pipettes can be used for many applications in the modern laboratory.

- Ergonomic design
- Colour coded
- Durable construction
- Thermally insulated
- Fully autoclavable
- User-friendly calibration
- Simple maintenance
- Made in Germany

SINGLE CHANNEL

SAVE OVER 50%!

WATCH NOW

PRODUCT CODE	VOLUME RANGE	PACK SIZE	PRICE	SALE
0020011210I	0.1 - 2µL	Each	LOGIN	LOGIN
0020011211I	0.5 - 10µL	Each	LOGIN	LOGIN
0020011213I	2 - 20µL	Each	LOGIN	LOGIN
0020011214I	10 - 100µL	Each	LOGIN	LOGIN
0020011215I	20 - 200µL	Each	LOGIN	LOGIN
0020011216I	100 - 1,000µL	Each	LOGIN	LOGIN
0020011217I	0.5 - 5mL	Each	LOGIN	LOGIN



SAVE
BIG!

Supplied with 3 ergonomic grips



soft grip

firm grip
large

firm grip
small

8 & 12 CHANNEL


WATCH NOW

PRODUCT CODE	VOLUME RANGE	PACK SIZE	PRICE	SALE
0020017943I	8 Channel 0.5 - 10µL	Each	LOGIN	LOGIN
0020017944I	8 Channel 10 - 100µL	Each	LOGIN	LOGIN
0020017945I	8 Channel 20 - 200µL	Each	LOGIN	LOGIN
0020017946I	12 Channel 0.5 - 10µL	Each	LOGIN	LOGIN
0020017947I	12 Channel 10 - 100µL	Each	LOGIN	LOGIN
0020017948I	12 Channel 20 - 200µL	Each	LOGIN	LOGIN



ACCESSORIES

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0020018987I	Pipette Rack Linea	Each	LOGIN	LOGIN
0020107134I	Pipette Stand Ronda	Each	LOGIN	LOGIN
0020016449I	Filter for Pette 0.5-5mL	Each	LOGIN	LOGIN
0020016450I	Filter for Pette 1-10mL	Each	LOGIN	LOGIN



0020018987I

0020107134I

Oven

125 BASIC

125L oven ideal for tempering, drying, aging and heating in the laboratory. Featuring an odourless insulation design, that also ensures lower operating costs.

- Fast heat up time
- Precise temperature control
- Clear display and control menu

Specifications:

- Temperature range: RT +5°C to 250°C
- 2 chrome-plated wire grid trays included, can fit up to 6 grid shelves
- Key lock function
- Glass door
- Overall dimensions: 700 x 825 x 650mm

PT100.31 temperature sensor (optional)

- Stainless steel
- Suitable up to 300°C



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0020003956I	Oven 125 Basic	Each	LOGIN	LOGIN
0020009564I	PT100.31 Temperature Sensor	Each	LOGIN	LOGIN

Incubator

INC 125 F DIGITAL

NEW PRODUCT

The INC 125 F digital convection incubator is ideal for microbiological applications up to 100°C with precise temperature control.

- Easy to clean
- Optimal temperature control
- Simple and intuitive operation
- Versatile accessories

Specifications:

- Temperature range: RT +8°C to 100°C
- 125L capacity
- Temperature stability: $\pm 0.3K$
- Dimensions: 700 x 720 x 635mm
- 2 chrome-plated wire grid trays included, can fit up to 6 grid shelves



IKA®

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0020117172I	Incubator INC 125 F Digital	Each	LOGIN	LOGIN

Homogeniser

T25 HOMOGENISER PACKAGE

High-performance dispersing instrument for volumes from 1 – 2,000 mL (H₂O) with digital speed display. It offers a wide speed range from 3,000 – 25,000rpm that enables users to work at high circumferential speeds even with small rotor diameters.

A broad choice of dispersing elements adds versatility. Applications range from homogenising wastewater samples, use in laboratory reactors, dispersion tasks under vacuum / pressure and sample preparation in medical research.

- Digital speed display
- Electronic speed control
- Electronic overload protection
- Stainless steel dispersing elements can be cleaned quickly and easily
- Plastic disposable dispersing elements are available in two sizes
- Error code display
- Quiet operation

Specifications:

- Volume range (H₂O): 1 - 2,000mL
- Max. viscosity: 5,000mPas
- Speed range: 3,000 - 25,000rpm
- Speed adjustment: Stepless
- Extension arm diameter: 13mm
- Extension arm length: 160mm
- Dimensions (WxHxD): 87 x 271 x 106mm
- Protection class: IP20
- Weight: 2.5kg

Supplied with:

- T25 Digital
- S25N/18G dispersing tool
- R182 bosshead clamp
- R1826 stand



Applications:


- Food & Beverage
- Pharmaceutical/Biotechnology
- Cosmetics & Personal Care
- Dairy
- Sample Preparation
- Scientific Research
- Environmental Analysis

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
00090242001	Homogeniser T25 Package	Each	LOGIN	LOGIN

IKA® Clearance Specials

Stirrer Magnetic I-MAG Package

- Powerful yet compact
- Stirs up to 300L
- Drive with very high IP protection class 64
- Electronically controlled motor
- Secure and ergonomic mounting



+ **SP300** frame included

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0020014226I-PACKAGE	Stirrer Magnetic I-MAG (with SP300 frame)	Each	LOGIN	LOGIN



Rotary Evaporator RV 8 V

- Precise manual lift handling
- Heating bath volume: 4L
- High-efficiency condenser with 1,600cm² cooling surface
- Integrated safety lift-out function
- Adjustable immersion angle

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0010003482I	Rotary Evaporator RV 8 V	Each	LOGIN	LOGIN



Tube Mill 100 Control

- Batch mill with disposable grinding chambers
- Adjustable timer from 5 seconds - 3 minutes
- USB interface
- Infinitely variable speed (5,000 - 25,000rpm)
- Maximum volume of 100mL

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0020007144I	Tube Mill 100 Control	Each	LOGIN	LOGIN



Tube Drive Homogeniser P Control Workstation

- Disposable tubes to prevent cross-contamination
- Programmable test conditions
- Adjustable change in direction of rotation
- Turbo button (speeds up to 8,000rpm)

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0025005836I	Tube Drive Homogeniser P Control Workstation	Each	LOGIN	LOGIN



Shaker Orbital MS3 Digital

- Compact, universal small shaker
- Wide range of attachments
- Attachment detection
- Timer with countdown function
- Continuous or touch operation (with standard attachment)

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0003319000I	Shaker Orbital MS3 Digital	Each	LOGIN	LOGIN



IKA® Clearance Specials

Stirrer Magnetic Hotplate RT10: Multi Position

- 10-position array
- Surface temperature infinitely adjustable up to 120°C
- Homogeneous temperature distribution
- Wear-free magnetic coils
- Speed range of 0 - 1,000rpm

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0003691100I	Stirrer Magnetic Hotplate RT10	Each	LOGIN	LOGIN



ThermoMix Matrix Orbital

- All-rounder, ideal for mixing, heating and cooling
- Large, clear display for intuitive operation
- Firm quick-release fastener
- Quick and effective mixing of liquids with up to 3,000rpm and 3mm shaker diameter
- Fast heating and cooling

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0030000731I	ThermoMix Matrix Orbital	Each	LOGIN	LOGIN



Stirrer Magnetic myPlate

- Lightweight and compact
- 135mm plate diameter
- Speed range of 100 - 2,500rpm
- Stirring volume of 1.5L
- Interchangeable and personal frame



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0025005156I	Stirrer Magnetic myPlate	Each	LOGIN	LOGIN



Overhead Stirrer Eurostar 40 Digital

- Digital speed display
- Up to 25L stirring volume (H₂O)
- Infinitely adjustable speed (0/30 - 2,000rpm)
- Push-through agitator shafts
- Overload protection
- Quiet operation

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
0004444000I	Overhead Stirrer Eurostar 40 Digital	Each	LOGIN	LOGIN

