

# Diaphragm Liquid Pump | Alligator 200




## FEATURES + BENEFITS

- **Time Saving**
- **Space Saving**
- **Efficient**



SPECIFICATIONS	ALLIGATOR 200
MAXIMUM VACUUM	313mBar
MAXIMUM FLOW RATE	45L/min
POWER INPUT	DC24V, 1.8A
MOTOR	60W Brushless
WEIGHT	2.4kg
WARRANTY	1 Year

See the Alligator 200 in action below:



PRODUCT CODE	PACK SIZE	PRICE	SALE
<a href="#">181200-02</a>	Each	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>

# Pipette Tips

## FEATURES + BENEFITS

- Made from 100% virgin polypropylene in fully automated class 100K clean room
- Swiss-made diamond finish mould for flawless, super smooth internal surface
- High transparency index
- Certified free from Dnase, Rnase, Endotoxin, Pyrogen & Human DNA & PCR Inhibitors
- Compatible with a wide range of pipettors
- Fully autoclavable



### Bulk Pack

Packaged loose in welded polyethylene bags

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	BREAK QTY
<a href="#">550.002.500</a>	0.1 - 10µL	PK1,000	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.510</a>	1 - 200µL	PK1,000	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.515</a>	100 - 1,000µL	PK1,000	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.527</a>	100 - 5,000µL	PK100	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.528</a>	500 - 10,000µL	PK100	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>

5,000ul & 10,000ul tips compatible with: Eppendorf, Gilson, IKA®, LabCo® & Rainin. Not suitable for Brand.



### Filtered Bulk Pack

- Specially developed PE filters stop aerosols & contamination from entering pipettor
- Packaged loose in welded polyethylene bags

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	BREAK QTY
<a href="#">550.002.600</a>	Filtered 0.1 - 10µL	PK1,000	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.605</a>	Filtered 1 - 20µL	PK1,000	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.610</a>	Filtered 1 - 200µL	PK1,000	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.615</a>	Filtered 100 - 1,000µL	PK1,000	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>



### Tip Box

- Autoclavable
- 8 x12 configuration
- Supplied empty

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	BREAK QTY
<a href="#">550.002.700</a>	Box for Pipette Tip 10µL	Each	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.710</a>	Box for Pipette Tip 200µL	Each	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.715</a>	Box for Pipette Tip 1,000µL	Each	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>



# Pipette Tips

## Racked, Non-Sterile Tips

- Tips are supplied pre-inserted into 96 place (8x12) racks
- Hinged boxes open to 120° for easy pipettor access


PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	BREAK QTY
<a href="#">550.002.535</a>	Racked 1 - 200µL	PK96 x 10	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.542</a>	Racked 100 - 1,000µL	PK96 x 6	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>



## Reloading Tips

- Reload rack with minimal effort
- Tips pre-arranged in 8x12 format
- Eco-friendly with minimal plastic consumption

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	BREAK QTY
<a href="#">550.002.570</a>	Reloading 0.1 - 10µL	PK96 x 10	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.575</a>	Reloading 1 - 200µL	PK96 x 10	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.580</a>	Reloading 1 - 1,000µL	PK96 x 5	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>



## Racked, Sterile Tips

- Tips are supplied pre-inserted into 96 place (8x12) racks
- Hinged boxes open to 120° for easy pipettor access
- Racked, sterile pipette tips are supplied pre-sterilised by irradiation

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	BREAK QTY
<a href="#">550.002.560</a>	Racked Sterile 1 - 200µL	PK96 x 10	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.567</a>	Racked Sterile 100 - 1,000µL	PK96 x 6	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>



## Low Retention, Filtered, Racked and Sterile Tips

- Specially developed PE filters stop aerosols & contamination from entering into pipettor
- Free from chemical & silicone additives = no leaching
- Minimum retention = maximum accuracy
- Manufactured from specially developed resins combined with high quality diamond finish mould creating super hydrophobic internal surface & ensuring maximum recovery of sample yield

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	BREAK QTY
<a href="#">550.002.670</a>	Low Retention Filtered Racked 10µL	PK96 x 10	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.675</a>	Low Retention Filtered Racked 200µL	PK96 x 10	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>
<a href="#">550.002.680</a>	Low Retention Filtered Racked 1,000µL	PK96 x 6	<a href="#">LOGIN</a>	<a href="#">LOGIN</a>

