

LabCo® Scientific Launch Specials



NITRILE GLOVES

MICROTOUCH® **Ansell**



- Blue Nitrile, Powder Free, Exam Glove
- Advanced stretch formulation for the ultimate comfort, fit and feel
- Tested for use with over 55 chemicals
- Ideal glove for almost all clinical settings and moderate risk examination procedures

PRODUCT CODE	SIZE	PACK SIZE
313041060	X SMALL	PK200
313041065	SMALL	PK200
313041070	MEDIUM	PK200
313041075	LARGE	PK200
313041080	X LARGE	PK200

*Order in multiples of 10 to receive full cartons

LATEX GLOVES

MICROTOUCH® **DermaClean® Ansell**



- Latex, Powder free, Exam Glove
- Secure grip thanks to its textured outer surface
- Easy donning with smooth inner surface
- Superior elasticity




PRODUCT CODE	SIZE	PACK SIZE
4570	X SMALL	PK100
4572	SMALL	PK100
4574	MEDIUM	PK100
4576	LARGE	PK100
4578	X LARGE	PK80

*Order in multiples of 10 to receive full cartons

labco THERMOMETER LABCO® DIGITAL MAX/MIN

- Thermometer & Humidity Meter
- Two probes to measure indoor/outdoor
- Memory Function for max/min date and time
- Calendar display
- Alarm clock function
- Dimensions 106 x 98 x 22mm
- 1 x AAA 1.5V Battery (supplied)



TECHNICAL SPECIFICATIONS	
TEMPERATURE RANGE X RESOLUTION	-50°C to 70°C x 0.1°C
TEMPERATURE ACCURACY	± 1°C
HUMIDITY RANGE X RESOLUTION	10% to 99% RH x 1% RH
HUMIDITY ACCURACY	± 5% RH

PRODUCT CODE	DESCRIPTION
750.500.050	Thermometer & Humidity Meter



labco ELECTRONIC TIMER UP/DOWN

- Free standing on flat surface
- Clips on to a pocket
- Magnetic back
- Easily readable digital display
- Programmable from 23 hours, 59 minutes, 59 seconds
- Two modes; Timer and Clock
- Count-up and count-down functions
- Memory function
- Alarm sounds for 1 minute
- 1 x AAA 1.5V Battery (supplied)



PRODUCT CODE	DESCRIPTION
LABTR01	Timer Electronic Up/Down

*1 FREE timer added automatically at checkout if eligible

Login or Register to view pricing
It's quick and easy!

- Check pricing & stock availability, anytime!
- Access to Online Exclusive Promotions
- Early Access to sales!
- Place orders
- Plus lots more!

labco FOUR CHANNEL TIMER

- Free standing on flat surface
- Clips on to a pocket
- Magnetic back
- Easily readable digital display
- Four independent countdown timer channels as well as a clock function
- Programmable from 99 hours, 59 minutes, 59 seconds
- Memory recall function
- Alarm sounds at 1 minute
- 1 x A76 Battery (supplied)



PRODUCT CODE	DESCRIPTION
PTR388	Four Channel

CONTAINS BUTTON BATTERY


This professional equipment is not intended for sale to the general public and/or to be used where children are present

LabCo pH INDICATOR STRIPS

Quick, easy and safe pH determination

Features

- Non Bleeding
- The indicator dyes in the test strips are chemically bound to the test pads. This patented technology safely prevents bleeding of dyes into your sample
- Long plastic handle
- Brilliant colour scale
- 6 x 85mm Strips
- Made in Germany



PRODUCT CODE	STR-PH-060	STR-PH-014
RANGE	pH 0 - 6	pH 0 - 14
GRADUATIONS	0, 0.5, 1.0, 1.5, 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14
PACK SIZE	PK100	PK100

TPS pH METER AQUA-pH

The Aqua-pH is your go-to choice for testing pH and Temperature, with the option for ORP as well.

The simple to use interface means it's easy to get testing, quickly and the fold out stand makes for hands free convenience.

The Aqua-Series are used from mines to meat testing, councils to construction and water treatment to wine making, plus more!

Features

- Sensors and buffers included
- Ready to use out of the box
- Waterproof to IP65
- 9V replaceable battery
- One touch calibration
- Fold out stand
- 1m Cable
- Optional waterproof carry case
- Australian Made



MODE	RANGES	RESOLUTION	ACCURACY
Temperature	-10.0 to 120.0°C	0.1°C	± 0.2°C
pH	0 to 14.00	0.01	± 0.01
ORP (mV)	0 to ± 1000.0	1	± -1

PRODUCT CODE	DESCRIPTION
121112-1TPS	Meter Aqua-pH pH/MV/ Temp 1 Meter Cable

Included in the kit:

- Aqua-pH meter
- pH Sensor
- Temp. sensor
- Buffer pH 4 & 7
- Quick Start Guide



labco FIXED HEIGHT TROLLEYS

Features

- Easy to assemble
- 100mm Casters (two with brakes)
- High strength plastic shelves 715 x 455mm
- Load capacity: 275kg

Login to view pricing now

PRODUCT CODE	DESCRIPTION	PACK SIZE
800.010.020	2 Shelves	EACH
800.010.030	3 Shelves	EACH



AND HR-250AZ COMPACT ANALYTICAL BALANCE

Compact and robust design. Fast and accurate precision measurement.

From 252g x 0.0001g or 3,889 grain x 0.002 grain

Features

- Fast stabilisation in 2 seconds
- Large removable breeze break with anti-static coating
- Compact foot print: 198 x 294mm
- Rotary sliding doors requiring no extra space at the rear
- Automatic calibration with an internal mass
- Clear, reverse-backlit LCD
- Fourteen weighing units including a user-programmable unit
- GLP/GMP/GCP/ISO compliant output
- Counting mode with the Automatic Counting Accuracy Improvement (ACAI) function
- Percent mode and statistical calculation function
- Adjustable response characteristics and hold function
- Comparator function with buzzer
- Auto power on/off function



PRODUCT CODE	DESCRIPTION	PAN SIZE
HR250AZ	252g X 0.1mg	90mm Ø

Supplied with large breeze break and AC adaptor
Industry Leading 5 Year Warranty

HUGE savings on balances!

AND FZ-i SERIES PRECISION BALANCES

The FZ-i Precision Balance brings top of the range performance to the price conscious. The easy, one touch internal calibration system provides high reliability and fail safe calibration.

Features

- One touch calibration
- Fast stabilisation time of 1 second
- Shock absorber function
- Complies with GLP, GMP, GCP and ISO

PRODUCT CODE	DESCRIPTION	PAN SIZE
FZ500I	520g x 0.001g	130mm Ø
FZ5000I	5200g x 0.01g	150mm Ø

Supplied with small breeze break, and AC adaptor and 3-year warranty.



AND EJ SERIES COMPACT SCALE

These accurate, high quality, attractive, stackable units are perfect for applications where portability or easy storage are paramount.

Features

- Stackable for easy storage
- Multiple measurement units including counting and percentage mode
- Security feature to protect against theft
- Optional density determination kit and breeze break
- Underhook weighing (for 0.1g model only)

PRODUCT CODE	DESCRIPTION	PAN SIZE
EJ610	610g x 0.01g	110mm Ø
EJ6100	6,100g x 0.1g	127 x 140mm Ø

Supplied with AC adaptor and 3-year warranty.



AND HT500 COMPACT SCALE

A low cost high quality solution to general weighing applications. Highly portable with long battery life (optional AC adaptor), the HT is perfect for field work applications.

Features

- Stackable case for easy storage
- Long battery life (approx. 450 hours - 4 x AA batteries, AC adaptor optional)
- Large LCD display
- Multiple weighing units (g, oz, ozt, lb-oz, tael & tola)
- Security slot for anti-theft lock

PRODUCT CODE	DESCRIPTION	PAN SIZE
HT500	510g x 0.1g	132 x 130mm

Supplied with perspex carry case and 1-year warranty.



LabCo® Scientific Launch Specials



labco ROTISSERIE TUBE MIXER

End over end, rotisserie action gives a gentle but effective mix

Features

- Speed range 10 to 70rpm
- Angle 0 to 90°
- Digital control with timer
- Wide range of attachments available (supplied with attachment for 16x50mL)



PRODUCT CODE

[400.104.040](#)

labco ORBITAL SHAKER 7.5KG

Compact shaker with orbital motion

Features

- Speed range 100 to 500rpm
- Orbital diameter 10mm
- LCD display for speed and timer
- Overspeed detection and protection
- Wide range of attachments available (supplied separately)



400.101.150 pictured with 400.101.157 attachment (not included in bonus offer)

PRODUCT CODE	DESCRIPTION
400.101.150	Orbital Shaker 7.5kg
400.101.155	Universal Attachment



400.101.150

**BONUS
ATTACHMENT**



400.101.155

LabCo[®] Scientific Launch Specials



labco MICRO CENTRIFUGE

Features

- 7,000rpm (2,680 x G)
- For 1.5mL, 2mL, 0.5mL, 0.2mL and PCR strip tubes

Supplied with:

- 400.003.010 Rotor (1.5/2mL x 8)
- 400.003.011 Rotor (0.2mL PCR tube strips x 2)



PRODUCT CODE	DESCRIPTION
400.003.005	Micro Centrifuge

labco MICRO LOW SPEED CENTRIFUGE

Ideal for routine medical practice in hospitals, clinics and research institutes.

Features

- Built-in fixed angle rotor
- Speed switch (rpm/rcf)
- LCD display
- Automatic lid-lock

TECHNICAL SPECIFICATIONS	
Maximum Speed	5,000rpm
Speed Range	300 to 5,000rpm
Maximum RCF	2,350 x G

PRODUCT CODE	DESCRIPTION
400.003.020	Centrifuge Micro Low Speed



labco CLINICAL CENTRIFUGE

Features

- Maximum speed: 4500rpm/ 2490 x G
- Rotor can be configured to spin multiple tube types (various adapters available)
- Door interlock
- Over speed and temperature detection
- Automatic internal diagnosis
- Short run time and Sound alert function
- Gentle braking at low speeds ensures efficient separation
- Parameters can be altered once desired speed is reached
- Short spin function, by holding the PULSE button
- Supplied with rotor for 15mL tubes (400.003.309)

PRODUCT CODE	DESCRIPTION	PACK SIZE
400.003.300	Clinical Centrifuge	Each

Supplied with rotor for 15mL tubes (400.003.309)



ROCKER 300 VACUUM PUMP

Login or Register to view pricing now

Features

- No air pollution, maintenance free
- Overflow protection
- Moisture trap with filter cartridge
- Vacuum regulator
- Thermal protection device



TECHNICAL SPECIFICATIONS	
VOLTAGE/FREQ.	220V / 50Hz
MAX. POWER (W)	60
MAX. CURRENT (A)	0.3
MAX. VACUUM (ULTIMATE)	105mbar / 26.8inHg
MAX. FLOW RATE	20L/min
MOISTURE TRAP/VACUUM REGULATOR	Yes
OVERFLOW PROTECTION	Yes
MOTOR ROTATION SPEED	1450rpm
HORSE POWER	1/8 HP
NOISE LEVEL	50dB
HOSE BARB	5/16" (8mm)
NET WEIGHT	4.1kg
DIMENSIONS (L X W X H)	26.8 x 13.5 x 20.4 cm
WARRANTY	2 years
PRODUCT CODE	167300-22

IKA® RCT BASIC HOTPLATE STIRRER

Features

- Integrated temperature control
- Exact temperature & speed setting via digital display, even when switched off
- Hot Top indicator >> hot surface warning to prevent burns!
- With adjustable safety circuit of heating plate temperature (50 to 360 °C)
- Bushing according to DIN 12878 for connecting a contact thermometer, e.g. ETS-D5, enables precise temperature control
- Enclosed assembly (IP42) guarantees long service life
- Includes PT1000 temperature sensor and protection cover



PRODUCT CODE	DESCRIPTION
00200026201	RCT Basic Stirrer

TECHNICAL SPECIFICATIONS	
STIRRING VOLUME (H ₂ O)	20L
SPEED RANGE	50 - 1,500rpm
PLATE DIMENSIONS	135mm Ø
PLATE MATERIAL	Aluminium Alloy
TEMPERATURE RANGE	Ambient - 310°C

Super Specials

Log in or Register for an online account to view prices!

lab_{co} MAGNETIC STIRRER 1.5L

A small and convenient, low foot print 1.5L magnetic stirrer to add to your laboratory!

- Square PET plate (120mm x 115mm)
- Stirring volume: 1.5L
- Speed range: 300 - 2,000rpm
- IP21 Protection
- Comes with electrode holder



PRODUCT CODE

[400.100.002](#)

COMPACT SCALE

The HT series brings optimum weighing performance and reliability, with an affordable price point.

- 510g x 0.1g
- Pan size: 132mm x 130mm
- Stackable
- Long battery life (~450 hours)
- Large LCD display
- Comes with a perspex carry case



PRODUCT CODE

[HT500](#)

lab_{co} CENTRIFUGE TUBES 15mL

Certified Non-Pyrogenic, DNase/RNase Free

- Sterile - Bulk Pack
- Conical Bottom
- 17mm (D) x 118mm (L)
- Autoclavable to 121°C
- Graduated from 1.5mL - 14.5mL x 0.5mL
- Max RCF 10,000g



PRODUCT CODE

[650.550.015](#)

PACK SIZE

PK500

lab_{co} pH ELECTRODE

Suitable for a variety of pH meters with a BNC connection. If you've been putting off getting a new pH electrode, now's your chance!

- Range: 0 - 14
- Temperature: 0 - 100°C
- Reference type: Ag/AgCl, Gel Filled, Single Junction
- Body type: Epoxy
- Connection: BNC



PRODUCT CODE

[950.051.010](#)

LABCO® PIPETTERS

One of our best sellers in the LabCo® range. Used by so many laboratories all around Australia!

- Fully autoclavable
- Lightweight & ergonomic
- Easy to calibrate
- Clear display

PRODUCT CODE	VOLUME RANGE
550.002.005	0.5 - 10µL
550.002.007	2 - 20µL
550.002.009	10 - 100µL
550.002.011	20 - 200µL
550.002.013	100 - 1,000µL
550.002.015	1,000 - 5,000µL

Login or Register to view prices



REFRACTOMETERS BRIX

Simplicity, efficiency and accuracy all in one handy device. Whether you're working with industrial fluids, oils, fruit juice, liquid sugars, beer and more, this refractometer is the perfect fit for your needs.

- ATC (Automatic Temperature Compensation)
- In-built calibration knob



PRODUCT CODE	APPLICATION	RANGE	RESOLUTION	ACCURACY	ATC
RHBN-10ATC	Brix, Industrial fluid	0-10% Brix	0.1%	±0.1%	10-30°C
RHBN-18ATC	Brix, Industrial fluid	0-18% Brix	0.1%	±0.1%	10-30°C
RHBN-32ATC	Brix, Industrial fluid	0-32% Brix	0.2%	±0.2%	10-30°C
RHBN-50ATC	Brix	0-50% Brix	0.5%	±0.5%	10-30°C
RHBN-62ATC	Brix	28-62% Brix	0.2%	±0.2%	10-30°C
RHBN-82ATC	Brix	45-82% Brix	0.5%	±0.5%	10-30°C

Login or Register to view prices